

EVOTIP

Loading protocol for
Evotip Pure




Evetips

Evetips are disposable, purification and loading trap columns.

Evetips fit most standard pipettes and commercial robotic systems.

Using Evetips results in simplified workflows that removes several sample handling steps and reduces injection cycle overheads.

Carry-over is minimized by eluting with a partial gradient in a very precise and reproducible volume that leaves impurities and high-molecular-weight contaminants behind.



Technical details

Specification	Value
Volume	0.5 - 200 μ l
Max. loading capacity	1 μ g of complex mixture or 1 pmol of simple protein digest

Ordering information

Part number	Value
EV2011	Evetip Pure, 1 box (96 tips)
EV2013	Evetip Pure, 10 boxes (960 tips)
EV2015	Evetip Pure, 50 boxes (4800 tips)
EV2018	Evetip Pure, 100 boxes (9600 tips)

Required Material

1. Box of Evotips.
2. Mass Spec Grade 1-Propanol, alternatively 2-Propanol.
3. Solvent A: Mass Spec Grade Water with 0.1% Formic Acid, JT Baker P/N: JT-9834-2 or equivalent.
4. Solvent B: Mass Spec Grade Acetonitrile with 0.1% Formic Acid, JT Baker P/N: JT-9832-2 or equivalent.
5. Multichannel pipette.
6. Solvent containers for multichannel pipette, e.g. Dual Solution Reservoir HEA20281A from Heathrow Scientific www.heathrowscientific.com/dual-solution-reservoir-i-hea20821a.
7. Single-well plate, e.g. 242811 from Thermo Scientific, www.thermofisher.com/order/catalog/product/242811. Do not use this plate with Solvent B (Acetonitrile).
8. Rack Holder w/bars (EV1166, Evosep).
9. Gloves.
10. Table centrifuge.
11. Fume hood.

Important info

- **At all centrifugation steps use max acceleration and de-acceleration.**
- **After each centrifugation step, empty liquid from the box.**
- **Please note that the following settings are guidelines, adjustment might be needed depending on the actual centrifuge.**
- **Centrifugation should be performed at room temperature.**



Follow local safety guidelines incl. the use of safety glasses, lab coat, gloves and a chemical safety cabinet when following this procedure.

Prepare solvents and containers

1. **Pour Solvent A and Solvent B into their respective solvent containers, and Propanol into the single-well plate.**

(The recommended dual purpose containers can either be used as a single reservoir or as 12 individual ones).




Rinse

1. Transfer 20 μ l Solvent B to all Evtips.
2. Place the Evtips in the centrifuge with appropriate counter balance and centrifuge at 800 g for 60 seconds.

IMPORTANT

EMPTY LIQUID FROM THE BOX AFTER EACH CENTRIFUGATION STEP.




Condition


1. Place the Evtip adapter rack on top of the single-well plate with Propanol.
2. Place the rack of Evtips onto this adapter so that the tips contact the Propanol.
3. Soak for a minimum of 10 seconds and then visually inspect that Evtips are pale white.

VERY IMPORTANT

**ENSURE THAT ALL EVOTIPS
ARE PALE WHITE AFTER SOAKING**




Let the Evtotips soak until all are pale white



 *Not pale white*


White vs pale white

**This step might take
up 10-20 seconds**



Equilibrate

1. With the Evtips still in the adapter rack, transfer 20 μ l Solvent A to each Evtip.



2. Place the Evtips in the centrifuge with appropriate counter balance and centrifuge at 800 g for 60 seconds.



IMPORTANT

EMPTY LIQUID FROM THE BOX AFTER EACH CENTRIFUGATION STEP.

Load sample

1. Place the Egotips in the original tray and transfer the sample (20 μ l) to each Egotip.

IMPORTANT


DO NOT PLACE THE EVOTIPS IN THE PROPANOL ADAPTOR AT THIS STAGE!

Solvent A can be added to the bottom of the Egotip box, to keep the Egotips wet during sample loading.

2. Place the Egotips in the centrifuge with counter balance and centrifuge at 800 g for 60 seconds.


IMPORTANT

EMPTY LIQUID FROM THE BOX AFTER EACH CENTRIFUGATION STEP.



Wash

1. Transfer 20 μ l Solvent A to each Evtip.



2. Place the Evtips in the centrifuge with appropriate counter balance and centrifuge at 800 g for 60 seconds.




IMPORTANT

EMPTY LIQUID FROM THE BOX AFTER
EACH CENTRIFUGATION STEP.

Preservation

1. Transfer 100 μ l Solvent A to each Evtip.
2. Place the Evtips in the centrifuge with appropriate counter balance and centrifuge at 800 g for 10 seconds to keep Evtips wet.



Storage recommendation

1. Add enough Solvent A to the bottom of the tray so the bottom of the Evtips are submerged and store at 4°C with lid on until analysis.



Quick guide

EvoTip Pure (EV2011)



RINSE

Wash dry EvoTips with 20 μ l Solvent B and centrifuge at 800 g for 60 s.



CONDITION

Soak in Propanol until all the EvoTips are pale white.



EQUILIBRATE

Equilibrate soaked EvoTips with 20 μ l Solvent A and centrifuge at 800 g for 60 s.



LOAD

Load samples on wet EvoTips (20 μ l in Solvent A) and centrifuge at 800 g for 60 s.



WASH

Wash EvoTips with 20 μ l Solvent A and centrifuge at 800 g for 60 s.



WET

Transfer 100 μ l Solvent A and centrifuge EvoTips at 800 g for 10 s to keep them wet.



EVOSEPP

NOTES

1. Ensure EvoTips are continuously wet after soaking in Propanol.
2. Transfer EvoTips to Solvent A tray in loading step to avoid that they dry out.
3. Prepare EvoTips just prior to MS analysis, but for longer storage add Solvent A in the bottom of the box and store refrigerated with lid on.
4. Empty liquid from the box after each centrifugation step.